

# **1.03.16 Lapping Machine for Thin Sections**





### **DESCRIPTION**

Machine for automatic thinning at a specified thickness of object glass blades and rocks glued to glass blades (thin blades).

It allows the heels to be used before sticking them to the glass blades.

# **PARTICULARITIES**

This machine is specialized in the automatic preparation of thin rock blades by rolling. It allows up to 12 samples to be taken simultaneously. The final thickness is adjusted using 3 micrometers graduated in 1/100 mm.

This package has unique advantages:

## Load the custom heels.

Heels can be very variable in thickness. They are always guided individually and the support load on the tray is individualized and applied in their center.

### Reliable and adjustable abrasive distribution in suspension.

The abrasive is kept permanently in suspension. The dispenser controlled by starting and stopping the tray allows the adjustment of the drip.

### Uniform distribution of the abrasive on the plate.


The abrasive is evenly distributed on the plate before coming into contact with the preparations.

### Permanent geometry: no disassembly of the tray is required.

The cast iron of the heavy-duty running-in plate guaranteed 20 years and its patented geometry holding device prevent disassembly for grinding throughout the life of the machine. **Patented sample holder.** 

A patented cleat allows the choice of the ground thicknesses.

## Easy maintenance and cleaning.



Updated in august 2020



BROTLAB.






# 1.03.16 LAPPING MACHINE FOR THIN SECTIONS

## STANDARD EQUIPMENT

- The machine with variable speed tray
- Patented cleat for finishing 12 blades 30x45 mm or 6 blades 45x60 mm (other dimensions on request)
- 1 Patented plate geometry support
- 1 Device with 3 or 6 cells for flat ground-in of the first face of the heels
- 2 Hands; 1 for the plate geometry retaining ring and abrasive splitter, 1 for the cleat and the heel
- 1 Built-in abrasive dispensing device
- 1 Water tube for holding samples in place by depression
  - ⇒ 1 Batch of start-up consumables:
  - $\Rightarrow$  1 Kg silicon carbide grain 800,
  - ⇒ 250 cm3 dispersing agent,
  - $\Rightarrow$  1 litre of passivating agent,
  - $\Rightarrow$  1 Pack of wipers,
  - $\Rightarrow$  2 Pissettes,
  - $\Rightarrow$  1 Brush for cleat,
  - $\Rightarrow$  1 Brush for the tray,
  - ⇒ 2 12-cell sample transport trays
- 1 Service Tool Lot



# **OPTION**

- Other voltage (to be specified at order)
- 1 Additional Hand (for second heel grounding device, allowing 2 sets of heels to be run in at the same time)

## **CHARACTERISTICS**

#### **Electrical**

Network: 230V single phase Output: 1 kW

### Capacity

Up to 12 thin sections 30x45 mm simultaneously.

# **Accessoiries**

- 1 Heel joint fitting (3 or 6) with ring+grid+weight (dimensions to be specified when ordering),
- 1 Vacuum pump with its water trap and its first maintenance batch,
- 1 Additional cleat,
- 1 1/100 analogue measuring stand
- 1 digital measuring stand at 1/100
- 1 digital measuring stand at 1/1000





